

*A Softer, Gentler Solution*

# Heelift® AFO



For three decades, the Heelift® name has been synonymous with quality pressure relieving products. Now, DM Systems, Inc., makers of Heelift Suspension Boot, introduces Heelift AFO, the newest member of the Heelift Family. Offering all of the benefits of a traditional ankle-foot orthosis, the new Heelift AFO also offers the quality and dependability one expects from a Heelift product.

## Indications

Heelift AFO is indicated for the following:

- Dorsiflexion ankle support for paralyzed legs
- Comfortable night splint
- Treatment/Prevention of plantar flexion contracture
- Treatment of plantar fasciitis
- Treatment/Prevention of foot drop deformity
- Treatment of Achilles tendonitis



Patent Pending

Authorized Code:  
L4396

**Heelift AFO offers a unique combination of benefits not found in other AFO products:**

- ▶ **Safety**
- ▶ **Customization**
- ▶ **Choice of Smooth or Convoluted Foam Interiors**
- ▶ **Ambulation with Assistance**
- ▶ **Washable/Autoclavable**





## Why Heelift® AFO?

### Safety

Unlike rigid orthotics that put patients at risk, the soft foam Heelift AFO prevents injury to the opposite leg.

The semi-rigid polypropylene brace is designed to provide the best possible control of the ankle's position and motion, while also allowing for limited weight bearing.

### Customization

Heelift AFO, like Heelift Suspension Boot, allows for patient-specific customization. The 5-centimeter, beveled, fixed pad elevates the heel in space and protects the Achilles tendon.

In addition, a spare pad with an adhesive backing is included to control external rotation.



**DM SYSTEMS, INC.**  
1316 Sherman Avenue  
Evanston, IL 60201 USA  
800.254.5438 • 847.328.9541  
info@dmsystems.com  
www.heelift.com

### Choice of Foam Interiors

Heelift AFO comes in a choice of smooth foam (ideal for edematous legs) or convoluted foam (a slightly cooler option).

### Ambulation with Assistance

Heelift AFO has a non-slip traction sole that allows limited weight bearing while preventing slippage. This feature is ideal for transfers from bed, wheelchairs, commode, or physical therapy. Assistance is recommended.

### Washable/Autoclavable

Like other products in the Heelift Family, Heelift AFO can be washed safely in a net or mesh laundry bag and can be autoclaved for sterility.



*There is Only One Heelift® AFO.*

To learn about Heelift AFO, visit  
[www.heelift.com/afo](http://www.heelift.com/afo).